


Salt Solution



Sterile Pyrogen Free Isotonic Balanced Salt Solution

Application

Irrigation of the anterior and posterior chamber during cataract surgery and other intraocular interventions, irrigation of alpha-chymotrypsin after enzymatic zonulolysis, intraoperative irrigation of the anterior chamber in cases of AC-bleedings. Irrigation of the posterior part of the eye during vitreoretinal surgery. Irrigation of the cornea during strabismus surgery and other extraocular surgeries, irrigation of the conjunctiva after fluorescein-instillation, irrigation of the eye and lacrimal ducts after etchings and burnings, moistening of the cornea and conjunctiva during laser treatments

Composition

sodium chloride, sodium acetate, sodium citrate, potassium chloride, calcium chloride, magnesium chloride, water for injection

Properties

The (chemical) composition, pH-value and osmolality of the BSS is due to the aqueous humor

Physico-chemical data

refractive index n_D^{20}	1.3344 - 1.3349
pH-value (25° C)	6.8 - 7.4
osmolality	300 ± 20 mOsmol

Packaging size

500 ml: Box with each 10 x 500 ml plastic bottles

Storage

at room temperature